

Gloves, Racks & General Lab Products

New Lower Prices on Extended Cuff Gloves

Nitrile 7 mil



Nitrile 4 mil



FREE with Purchase*

Box of Truffles or 3 Pack Chocolates with any order of \$200 - \$599

Flash Drive & Stylus Pen with any order of \$600 - \$999

Cooler Bag & Collapsible Bowl with any order of \$1000 & Over



PCR Reagents & Plastics FLYER 2018 OFFERS VALID UNTIL Dec. 28, 2018
Promo Code Q34-2018 Required at Time of Order

Promotional Giveaways will be withheld from certain companies/organizations due to their policies.
Diamed reserves the right to withhold these promotional items if the orders are not going direct to the end user.

*All promo free items available only with purchase
of consumable products listed in Flyer.

**Buy 9 Get ONE FREE
on all NITRILE Gloves**

*Mix'n Match
Sizes*
for Volume Discount



Nitrile 3.9 mil

- Silky Soft
- Finger Textured
- Powder Free
- 1.5 AQL
- 200 per box

**Best Value
200/Box**

*Mix'n Match
Sizes*
for Volume Discount

\$16.49 ea.
200/Box

Extra Small	Small	Medium	Large
DiaSILK-XSPF	DiaSILK-SPF	DiaSILK-MPF	DiaSILK-LPF



Nitrile 4.3 mil

- Superior Strength
- Fully Textured
- Powder Free
- 1.5 AQL
- 100 per box

*Mix'n Match
Sizes*
for Volume Discount

\$11.29 ea.
100/Box

Extra Small	Small	Medium	Large	Extra Large
TECNITE-XSPF	TECNITE-SPF	TECNITE-MPF	TECNITE-LPF	TECNITE-XLPP



Nitrile 3 mil

- Polymer Coated
- Finger Textured
- Powder Free
- 1.0 AQL
- 100 per box

CHEMO Tested

*Mix'n Match
Sizes*
for Volume Discount

\$14.99 ea.
100 / Box

Extra Small	Small	Medium	Large
Dia3009-XSPF	Dia3009-SPF	Dia3009-MPF	Dia3009-LPF



Nitrile 4 mil

- Comfortable Grip
- Fully Textured
- Powder Free
- 1.0 AQL Tear Proof
- 100 per box

*Mix'n Match
Sizes*
for Volume Discount

\$14.99 ea.
100 / Box

Extra Small	Small	Medium	Large
Dia4009-XSPF	Dia4009-SPF	Dia4009-MPF	Dia4009-LPF



Nitrile 5 mil

- Black Matte
- Finger Textured
- Beaded Cuff
- Powder Free
- 1.5 AQL
- 100/box

\$13.99 ea.
100 / Box



Nitrile 7 mil

- Extended cuff
- Finger Textured
- Powder Free
- 1.5 AQL
- 50/box

New Lower Price

*Mix'n Match
Sizes*
for Volume Discount

\$11.99 ea.
50/Box



Nitrile 4 mil

- Extended Cuff
- Finger Textured
- 12 Inches Long
- Powder Free
- 1.5 AQL
- 100/box

New Lower Price

*Mix'n Match
Sizes*
for Volume Discount

\$10.99 ea.
100/Box

Small	Medium	Large	Extra Large	Extra Small	Small	Medium	Large	Extra Small	Small	Medium	Large
DiaSTRONG-SPF	DiaSTRONG-MPF	DiaSTRONG-LPF	DiaSTRONG-XLPP	DiaSAFE-XSPF	DiaSAFE-SPF	DiaSAFE-MPF	DiaSAFE-LPF	DiaSOFTXC-XSPF	DiaSOFTXC-SPF	DiaSOFTXC-MPF	DiaSOFTXC-LPF

Need Speciality Gloves?

We have them!
Call for Details



GLOVES, RACKS, BOXES & MISCELLANEOUS



Latex 6 mil

- Polymer Coated
- Superior Strength
- Honeycomb Grip
- Beaded Cuff
- Powder Free
- 1.5 AQL
- 100/box

1 - 9 Boxes
\$14.99 ea.

Mix'n Match
Sizes

Buy 10 + for
\$13.99 ea.
100/Box

Extra Small	Small	Medium	Large	Extra Large
PRELG-XSPF	PRELG-SPF	PRELG-MPF	PRELG-LPF	PRELG-XLPF



Latex 5 mil

- Polymer Coated
- Fully Textured
- Beaded Cuff
- Powder Free
- 1.5 AQL
- 100/box

1 - 9 Boxes
\$9.99 ea.

Mix'n Match
Sizes

Buy 10 + for
\$8.99 ea.
100/Box

Extra Small	Small	Medium	Large	Extra Large
TECLG-XSPF	TECLG-SPF	TECLG-MPF	TECLG-LPF	TECLG-XLPF

Cotton Glove Liners

- Avoid Allergy Contamination

Small / Medium

TECCGL-1	12 Pair	\$40.30
TECCGL-2	144 Pair	\$288.29

Large

TECCGL-3	12 Pair	\$42.89
TECCGL-4	144 Pair	\$308.84



Biodegradable DiaECO

DiaECO Nitrile Glove attracts microbial activity in a landfill, promoting quicker decomposition

Buy 9 Get ONE FREE
on all NITRILE Gloves

Nitrile 4 mil

- 100% Nitrile
- 4 mil (0.10mm) thick
- 9.5 inches (240mm) long
- Fingertip Texturing
- AQL 1.5

Mix'n Match
Sizes

\$17.49 ea.
10/Box

Extra Small	Small	Medium	Large	Extra Large
DiaECO-XSPF	DiaECO-SPF	DiaECO-MPF	DiaECO-LPF	DiaECO-XLPF

DiaECO biodegrades in biologically-active landfills thus providing an effective single-use hand protection solution. The unique biodegradable characteristic is made possible by an innovative, patent-pending technology. This organic additive is designed to make DiaECO attractive to microbial activity in a landfill and upon its consumption, micro-organisms excrete enzymes that depolymerise the nitrile. The result is biogases and inert humus. Normal nitrile gloves take decades, if not hundreds of years, to biodegrade and break down in landfill whereas the biodegradation rate of DiaECO is far more rapid.

Testing by an independent research lab has confirmed that 20 per cent of a DiaECO glove will have biodegraded in 160 days. Importantly there is no difference in shelf life as DiaECO does not start to biodegrade until it is in landfill.

-80°C Safe Freezer Labels

MML-? 500 Sets of Lid & Tube Labels

\$83.95 /Set

Specify Colour:

B - Blue, G - Green, R - Red,
W - White & Y - Yellow

Mix'n Match Colours



Laboratory Tape - 500 Inches per Roll

Mix'n Match
Colours

		Price	10 Price
ST12-E?	1/2" Wide	\$3.09 ea.	\$26.95
ST34-E?	3/4" Wide	\$4.19 ea.	\$39.95
ST10-E?	1" Wide	\$5.19 ea.	\$44.95

Specify Colour:

B - Blue, G - Green, O - Orange, P - Pink,
R - Red, V - Violet, W - White, Y - Yellow



Cryo Dots

- 9.5mm Dots
- 5000 Dots
- Blue, Green, Purple, Red & Yellow (1000 ea.)
- Designed to withstand extreme cold -196° C

STDOTS-A5 **\$57.95**




Autoclave Tape

60 Yards per Roll

Lead-Free

ST010-EAT	3/4" Wide	\$22.49 ea.
ST020-EAT	1" Wide	\$30.95 ea.

**Special
Buy 2 Rolls
1" Wide
Get 1 FREE**



GLOVES, RACKS, BOXES & MISCELLANEOUS

Weighing Dishes

- White polystyrene
- Round balanced footprint
- Contoured Sides for Easy Secure handling
- Flexible for easy sample pouring
- Conventional Heating to 90 °C.

Cat. #	Shape	Size (mm)	Qty	Price
TECWB-1	Square with Round base	42 x 8 Deep	500	\$40.40
TECWB-2	Square with Round base	88" x 25 Deep	500	\$68.97
TECWB-3	Square with Round base	140 X 22 Deep	500	\$147.82



Antistatic Weighing Dishes

- Made of antistatic white polystyrene
- Will withstand temperatures up to 80 °C.

Cat. #	Shape	Capacity	Size (mm)	Qty	Price
DLAE500-1	Square	10ml	40 x 40 x 8 H	500	\$20.78
DLAE500-2	Square	100ml	78 x 78 x 25 H	500	\$48.25
DLAE500-3	Square	200ml	127 x 127 x 25 H	500	\$97.82
DLAE502-1	With Spout	12ml	50 x 37 x 8 H	250	\$21.69
DLAE502-2	With Spout	140ml	128 x 76 x 25 H	250	\$44.98
DLAE502-3	With Spout	270ml	180 x 117 x 25 H	250	\$70.32
DLAE504-1	Hexagon	9ml	25 Base, 38 Top x 10 H	500	\$21.98
DLAE504-2	Hexagon	58ml	47 Base, 70 Top x 20 H	250	\$44.48
DLAE504-3	Hexagon	203ml	85 Base, 115 Top x 22 H	250	\$75.55
DLAE504-4	Hexagon	355ml	95 Base, 130 Top x 30 H	250	\$89.61



Antistatic Weighing Dishes

- Made of antistatic polystyrene
- Will withstand temperatures -20 °C to 90 °C

Cat. #	Shape	Capacity	Size (mm)	Qty	Price
MTCB6001B	Square (Black)	20ml	45 x 45	500	\$29.05
MTCB6001W	Square (White)	20ml	45 x 45	500	\$31.47
MTCB6002B	Square (Black)	100ml	80 x 80	500	\$45.38
MTCB6002W	Square (White)	100ml	80 x 80	500	\$49.17
MTCB6003B	Square (Black)	250ml	140 x 140	500	\$88.95
MTCB6003W	Square (White)	250ml	140 x 140	500	\$96.37
MTCB6501B	Diamond (Black)	5ml	35 x 55	500	\$26.55
MTCB6501W	Diamond (White)	5ml	35 x 55	500	\$26.55
MTCB6502B	Diamond (Black)	30ml	65 x 85	500	\$42.87
MTCB6502W	Diamond (White)	30ml	65 x 85	500	\$42.87
MTCB6503B	Diamond (Black)	100ml	100 x 125	500	\$46.81
MTCB6503W	Diamond (White)	100ml	100 x 125	500	\$46.81



Disposable Funnel

- All sizes are tapered at a 60° angle with inside fluting

DLAG980-1	57 mm Top ID, 69mm High, Polystyrene	100	\$43.43
DLAG980-2	64 mm Top ID, 76mm High, Polystyrene	100	
DLAG980-3	56 mm Top ID, 69mm High, Polypropylene	100	
DLAG980-4	63 mm Top ID, 75mm High, Polypropylene	100	

Autoclave Bags

- With BioHazard Sign & Writing Area,
- Maximum Temperature 136° C
- Polypropylene
- Bulk 200/Pk.


MTCA9000C	12"x26" - Clear	\$137.67
MTCA9000R	12"x26" - Red	\$137.67
MTCA9002C	24"x32" - Clear	\$274.95
MTCA9002R	24"x32" - Red	\$274.95



Bench Liner

- Absorbent top layer with waterproof bottom layer
- 18 inch wide roll with perforations

AD500-BL1813	18"x13.5" - 40 Sheets/Roll	\$9.59
AD500-BL1827	18"x27" - 20 Sheets/Roll	\$9.59



GLOVES, RACKS, BOXES & MISCELLANEOUS

Convenience Items

Plastic Wrap

TECWRAP228	All Purpose Plastic Wrap, 30.48cm x 228.6M	1 roll	\$17.95
TECWRAP914	All Purpose Plastic Wrap, 30.5cm x 914.4M (comes in easy to use box with cutter)	1 roll	\$41.95

Aluminum Foil

TECFOIL30	Alcan Aluminum Foil, 30.5cm x 30.48M	1 roll	\$8.95
TECFOIL200	Alcan Aluminum Foil Plus , 30cm x 200M	1 roll	\$52.99
TECFOILHD100	Heavy Duty , Alcan Aluminum Foil Plus, 45cm x 100M	1 roll	\$49.95

Parchment Paper

TECNPAPER50	Non-Stick Parchment Paper, 216C/420F Incubator safe, 38cm x 49.9M	1 roll	\$13.62
-------------	--	--------	---------

Miscellaneous

TECSEAL22	Glad - Press'n Seal, 30cm x 21.6M (6.5m ²)	1 roll	\$13.49
TECSEAL43	Glad - Press'n Seal, 30cm x 43.4M (13m ²)	1 roll	\$28.95
TECTOWEL	Single Fold Paper Towels (250/Bundle)	3 Packs	\$16.95
TECSOAP	Soft Soap Refill - Soothing Aloe Vera	2.36 litre	\$13.95
TECMARKER	Sharpie Markers - Fine Point	5 Pack	\$13.95
TECZIP612	Reclosable Bags, Heavy Duty 4mil, 6x12"	100 Pack	\$32.95



Chromatography Columns

5 & 7" Columns

- Bulk, 200/Pk.
- Polyethylene Caps
- Include reservoir caps, end caps, and Polyethylene (PE) filter matrix supports

E/S208-3366-050	5" PS with coarse (90-130 um) filters	\$283.31
E/S208-3367-050	5" PS with medium 45-90 um) filters	\$283.31
E/S208-3383-060	5" PP with coarse (90-130 um) filters	\$270.52
E/S208-3384-060	5" PP with medium (45-90 um) filters	\$270.52
E/S208-3385-060	5" PP with fine (15-45 um) filters	\$270.52
E/S208-3393-060	7" PP with medium (45-90 um) filters	\$312.64

Universal Columns, 5.5" Long

- Bulk, 200/Pk.
- Polypropylene
- Chemically resistant to weak acid and alkali
- Total capacity of 25 mL

E/S208-3388-060	5½" PP with medium (45-90 um) filters	\$254.26
E/S208-3389-060	5½" PP with fine (15-45 um) filters	\$254.26




Cuvettes

- 100 per Tray, 1000/Pk.
- 10mm Path Length
- Optically Clear Polystyrene
- 45mm High
- Range: 340-800 nm

STK2400	3.5 ml Standard Cuvette, 2 Sides Clear	\$143.49
STK2405	3.5 ml Fluorimeter Cuvette, 4 Sides Clear	\$159.61
STK2410	1.5 ml Semi-Micro Cuvette, 2 Sides Clear 4mm Width	\$149.95
STK2401	Cuvette Caps, 2 x 500/Bag	\$45.47

- 100 per Tray, 500/Pk.
- 10mm Path Length
- UVT Acrylic
- Packaged in special trays that only contact the corners and protect the optical surfaces
- Ribbed cuvettes provide added stability of cuvette orientation and are recommended for kinetic studies

E/S211-3113-0A0	Semi-micro, 0.6 ml min volume, special ribbed design for added stability	\$218.75
E/S211-3127-0A0	Standard, 4 ml capacity, frosted two sides	\$206.08
E/S211-3128-0A0	Standard, 4 ml capacity, clear four sides	\$223.91
E/S211-3169-0A0	Semi-micro, 1.0 ml min volume, clear four sides, jacketed	\$231.00



GLOVES, RACKS, BOXES & MISCELLANEOUS

Call for Volume Pricing for Boxes & Dividers

Cardboard Storage Boxes & Dividers

- 5.25 x 5.25" Square, Large Box 5.75 x 5.75"

Boxes 2"

TECCB2	2" Box Only with Lid	ea.	\$2.95
TECCB2-81	2" Box with Lid & 81 Cell Divider	ea.	\$3.49
TECCB2-100	2" Box with Lid & 100 Cell Divider	ea.	\$3.95

Boxes 3"

TECCB3	3" Box Only with Lid	ea.	\$2.95
TECCB3-81	3" Box with Lid & 81 Cell Divider	ea.	\$3.95
TECCB3-100	3" Box with Lid & 100 Cell Divider	ea.	\$5.29
TECCB3-L	3" Large Box Only with Lid	ea.	\$2.95




Dividers

TECCD16	16 Cell Divider	\$2.75
TECCD16-L	16 Cell Divider, for CB3-L Box	\$5.50
TECCD36	36 Cell Divider	\$2.75
TECCD36-L	36 Cell Divider, for CB3-L Box	\$5.50
TECCD49	49 Cell Divider	\$2.75
TECCD64	64 Cell Divider	\$2.75
TECCD81	81 Cell Divider	\$2.29
TECCD100	100 Cell Divider	\$2.75

Speciality Box for Short Vials

E/S240-8225-000	Cardboard Freezer Box with 25 cell partitions for 2.5 ml large I.D. vials 5.25 x 5.25" x 1.25" H	12/Pk.	\$74.37
-----------------	--	--------	---------



Custom Boxes & Dividers Available

Microtube Freezer Racks

Lift Off Lid
81 Place



TECFR81-? \$9.49 ea.

**Buy 15
Get TWO
FREE***
See Page 8
for Details

Hinged Lid
100 Place



Not Available in Black
TECTR100-? \$9.95 ea.

Specify Colour: B - Blue, BK - Black, G - Green, N - Natural, O - Orange, P - Pink & Y - Yellow

15 & 50ml Tube Rack

- Operating Temperature: -90°C to +121°C (autoclavable)
- Unique pull tab
- Clear printed grid
- Fits in standard freezer racking
- Polypropylene

TECTR15-36G	Holds 36, 15ml Tubes	ea.	\$27.34
TECTR50-16B	Holds 16, 50ml Tubes	ea.	\$27.34



15ml & 50ml Combo Tube Rack

Holds 50 tubes

MTCR1070 \$75.90 (5/pk)



Test Tube Rack

- Lightweight, stackable and autoclavable
- All racks measure 246 mm x 104 mm x 64 mm

Blue	Red	White	Capacity	
E/S240-3200-B8SI	E/S240-3200-R8SI	E/S240-3200-W8SI	96 Place, Fits 13mm Tubes	\$27.75
E/S240-3201-B8SI	E/S240-3201-R8SI	E/S240-3201-W8SI	60 Place, Fits 17mm Tubes	\$27.75
E/S240-3204-B8SI	E/S240-3204-R8SI	E/S240-3204-W8SI	21 Place, Fits 30mm Tubes	\$27.75



30 & 50ml Tube Racks

E/S240-3500-G10	Polystyrene, Green	\$47.01
E/S240-3500-Y10	Polystyrene, Yellow	\$47.01

GLOVES, RACKS, BOXES & MISCELLANEOUS

Tube Racks - 80 Place

Accommodates all types of screw cap microtubes from 0.5 to 2 ml as well as 1 to 5 ml cryogenic vials. Flip it over and you can store up to 60 PCR or microcentrifuge tubes from 0.2 to 0.5 ml

- Size: 223 x 67 x 27 mm H (9 3/16 x 2 5/8 x 1 1/16 in. H)

DLAT1000-80?	80 Place Tube Rack	10/pk.	\$46.24
DLAT1002-80	Transparent Lid for Above	10/pk.	\$32.27

Specify Colour: B – Blue, G – Green, P – Pink, R – Red, Y – Yellow

Smaller model holds 25 tubes, or 16 PCR or microcentrifuge tubes

DLAT1000-25?	25 Place Tube Rack	10/pk.	\$34.24
--------------	--------------------	--------	----------------

Specify Colour: B – Blue, R – Red, Y – Yellow



**Buy 15
Get TWO
FREE***
See Page 8
for Details

GRIPPER Rack - Ideal for Waterbaths

This GRIPPER 50-place rack holds 12 x 75 mm and 13 x 100 mm tubes securely in place thanks to silicone tabs surrounding the base of each tube while in the rack.

- Size: 250 x 128 x 50 mm H (9 7/8 x 5 x 2 in. H).

DLAT1020-500	50 Place GRIPPER Rack	2/pk.	\$56.86
--------------	-----------------------	-------	----------------



Multi Purpose Racks

- Submersible & Stable
- Great for Waterbaths

Submersible

Half Size

Standard Size

DLAT1200-13?	84 Place	\$22.35 ea.	Up to 13mm Tubes	DLAT1220-13?	42 Place	\$13.02 ea.
DLAT1200-16?	60 Place	\$22.35 ea.	Up to 16mm Tubes	DLAT1220-16?	30 Place	\$13.02 ea.
DLAT1200-30?	18 Place	\$22.35 ea.	Up to 30mm Tubes	DLAT1220-30?	9 Place	\$13.02 ea.

Specify Colour: O-Orange, Y-Yellow, B-Blue, G-Green or L-Lilac



Freezer Racks - 50 Place

- Durable and robust Polypropylene
- Hinged Lid
- Stores tubes up to 2.0 ml.
- Suitable for temperatures as low as -90°C.

TECTR50-?	Hinged Freezer Rack	\$7.83
-----------	---------------------	---------------

Specify Colour:

B – Blue, O – Orange, Y - Yellow, N -Natural



PCR Tube Racks

Capacity:
96 x 0.2ml Tubes or
8 & 12 Strip Tubes

TECMR0296	\$9.39
-----------	---------------

Colour Choice Subject to Stock Availability




PCR Workup Racks & Workstation Assembly

Capacity:
24 x 1.5/2.0ml Tubes. 16 x 0.5ml Tubes &
96 x 0.2ml Tubes or Strips

TECMR241696	\$30.87 ea.
-------------	--------------------

Colour Choice Subject to Stock Availability



Includes Lid

64 Place Flipper Racks

Available Colours:
N – Natural, B – Blue BK - Black, G – Green,
R - Red, & V – Violet, (Orange Discontinued)

TECFR64-?	\$10.99 ea.
-----------	--------------------

Colour Choice Subject to Stock Availability



96 Place Flipper Racks

- 96 Place, 8 rows x 12 columns with indexing
- Polypropylene
- Operating Range: -90 to 121°C
- Clear Lids

TECMR96-?	\$14.66 ea.
-----------	--------------------

Specify Colour

Standard Colours:

B – Blue, BK – Black, G – Green, N – Natural,
O – Orange, R – Red, V – Violet, Y – Yellow

Fluorescent Colours:

FB – Blue, FG – Green, FO – Orange,
FP - Pink, FY – Yellow



Combination Flipper Racks

Capacity:
12 x 1.5/2.0ml Tubes & 96 x 0.2ml
Tubes or Strips or 96 x 0.5ml Tubes

TECMR12-96	\$51.90 ea.
------------	--------------------

Colour Choice Subject to Stock Availability



Lid Included



Microtube Racks - 80 Place

Specify Colours:
 B – Blue, BK – Black,
 G – Green, N – Natural,
 O – Orange, P – Pink,
 V – Violet & Y – Yellow

TECMR80-? \$4.59 ea.

Specify Colour



Container Only for 80 Place Rack

Natural Container
 Lid & Base Only
 Holds one TECMR80

TECTR80 \$12.83 ea.



Shown with
 TECMR80 Rack

4 Way Flipper Racks

Available Colours:
 N – Natural,
 Fluorescent Colours:
 FB – Blue FG – Green,
 FO – Orange, FP – Pink,
 & FY – Yellow

TECFR4-?? \$9.49 ea.



**Buy 15
 Get 2
 FREE***
**on All Racks
 Pages 6-8**
 Lowest Priced Racks are FREE*

4 Way Flipper Racks - with Tray

Tube Size Options:
 • (12) 0.5 ml tubes,
 • (9) 1.5-2.0 ml tubes,
 • (6) 15 ml tubes or
 • (2) 50 ml tubes

TECFR44 \$67.95 Set



1.5ml Microtube Rack

- 100 Place
- Lid Included

NLR8000 \$9.95 ea.



32 Place Flipper Racks

Fluorescent Colours:
 FB – Blue, FG – Green, FO – Orange,
 FP – Pink, FY – Yellow

TECFR32-8 \$13.39 ea.

Colour Choice Subject
 to Stock Availability



Timers

- Batteries Included



#DiaTIMER-VMCCC



#DiaTIMER-COCC

DIATIMER-MCCC	Traceable® Timer with Magnetic Clip	\$39.95
DIATIMER-COCC	Traceable® Clip On Timer	\$50.95



#DiaTIMER-D100CC

DIATIMER-VMCCC	Traceable® Vibrating Timer with Magnetic Clip/Stand	\$42.95	DIATIMER-D100CC	Traceable® 100 Hour Dual Timer	\$57.95
----------------	---	---------	-----------------	--------------------------------	---------



#DiaTIMER3CC



#DiaTIMER4CC

TRACEABLE® Timers
 Comply to ISO 1750.01
 & are Traceable to NIST
 Standards



#DiaTIMER



#DiaTIMER-LD

DIATIMER3CC	Traceable® 3-Channel Alarm Timer	\$75.95	DIATIMER	Laboratory Timer/Clock	\$28.42
DIATIMER4CC	Traceable® 4-Channel Alarm Timer	\$65.95	DIATIMER-LD	Large Display with Stainless Housing	\$28.42

You must use **Promotional Code Q34-2018**
 on **ALL** Purchase Orders to Receive
Special Flyer Pricing - No Exceptions
OFFERS VALID UNTIL December. 28, 2018

TOLL FREE: 1.800.434.2633
FAX: 1.800.293.2987
LOCAL: 905.625.6021
FAX: 905.625.6280
WEBSITE: www.diamed.ca
EMAIL: diamed@diamed.ca

FREE Shipping
 On Consumable Orders
 Over \$200.⁰⁰
 (Ground Only)



© Diamed Lab Supplies Inc. All rights reserved. Diamed, DiaTec, DiaSoft, DiaStretch, DiaGrip, TeclG, Tec Choice are trademarks or registered trademarks of Diamed Lab Supplies Inc. **Terms & Conditions:** Prices and promotions may change without notice. Taxes and shipping charges, if applicable, are not shown. Diamed Lab Supplies Inc. is not responsible for pricing or other errors and reserves the right to modify or cancel orders if such errors occur. **Promotional & Free Items:** Must be requested at time of order to receive. Void where prohibited by institutional policy. Qualifying amount calculated before freight and applicable taxes. **Payment Terms:** Institutional Purchase Order, Visa, MasterCard, and American Express are accepted.